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CellTiter-Glo[®] Luminescent Cell Viability Assay

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1. Crouch, S.P. *et al.* (1993) The use of ATP bioluminescence as a measure of cell proliferation and cytotoxicity. *J. Immunol. Methods* **160**, 81–8.
2. Farfan, A. *et al.* (2004) Multiplexing homogeneous cell-based assays. *Cell Notes* **10**, 2–5.
3. Riss, T., Moravec, R. and Niles, A. (2005) Selecting cell-based assays for drug discovery screening. *Cell Notes* **13**, 16–21.
4. Kangas, L., Grönroos, M. and Nieminen, A.L. (1984) Bioluminescence of cellular ATP: A new method for evaluating cytotoxic agents in vitro. *Med. Biol.* **62**, 338–43.
5. Lundin, A. *et al.* (1986) Estimation of biomass in growing cell lines by adenosine triphosphate assay. *Methods Enzymol.* **133**, 27–42.
6. Sevin, B.U. *et al.* (1988) Application of an ATP-bioluminescence assay in human tumor chemosensitivity testing. *Gynecol. Oncol.* **31**, 191–204.
7. Gerhardt, R.T. *et al.* (1991) Characterization of in vitro chemosensitivity of perioperative human ovarian malignancies by adenosine triphosphate chemosensitivity assay. *Am. J. Obstet. Gynecol.* **165**, 245–55.
8. Petty, R.D. *et al.* (1995) Comparison of MTT and ATP-based assays for the measurement of viable cell number. *J. Biolumin. Chemilumin.* **10**, 29–34.
9. Cree, I.A. *et al.* (1995) Methotrexate chemosensitivity by ATP luminescence in human leukemia cell lines and in breast cancer primary cultures: Comparison of the TCA-100 assay with a clonogenic assay. *AntiCancer Drugs* **6**, 398–404.
10. Maehara, Y. *et al.* (1987) The ATP assay is more sensitive than the succinate dehydrogenase inhibition test for predicting cell viability. *Eur. J. Cancer Clin. Oncol.* **23**, 273–6.
11. Stanley, P.E. (1986) Extraction of adenosine triphosphate from microbial and somatic cells. *Methods Enzymol.* **133**, 14–22.
12. Beckers, B. *et al.* (1986) Application of intracellular ATP determination in lymphocytes for HLA-typing. *J. Biolumin. Chemilumin.* **1**, 47–51.
13. Andreotti, P.E. *et al.* (1995) Chemosensitivity testing of human tumors using a microplate adenosine triphosphate luminescence assay: Clinical correlation for cisplatin resistance of ovarian carcinoma. *Cancer Res.* **55**, 5276–82

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产品	规格	目录号
ApoLive-Glo™ Multiplex Assay	10ml	G6410
ApoTox-Glo™ Triplex Assay	10ml	G6320
CellTiter-Fluor™ Cell Viability Assay (fluorescent)	10ml	G6080
CellTiter-Blue® Cell Viability Assay (resazurin)	20ml	G8080
CellTiter 96® AQueous One Solution Cell Proliferation Assay (MTS, colorimetric)	200 assays	G3582
CellTiter 96® AQueous Non-Radioactive Cell Proliferation Assay (MTS, colorimetric)	1,000 assays	G5421
CellTiter 96® AQueous MTS Reagent Powder	1g	G1111
CellTiter 96® Non-Radioactive Cell Proliferation Assay (MTT, colorimetric)	1,000 assays	G4000

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产品	规格	目录号
CytoTox-Glo™ Cytotoxicity Assay (luminescent)*	10ml	G9290
Mitochondrial ToxGlo™ Assay*	10ml	G8000
MultiTox-Glo Multiplex Cytotoxicity Assay (luminescent, fluorescent)*	10ml	G9270
MultiTox-Fluor Multiplex Cytotoxicity Assay (fluorescent)*	10ml	G9200
CytoTox-Fluor™ Cytotoxicity Assay (fluorescent)*	10ml	G9260
CytoTox-ONE™ Homogeneous Membrane Integrity Assay (LDH, fluorometric)*	200–800 assays	G7890
CytoTox-ONE™ Homogeneous Membrane Integrity Assay, HTP	1,000–4,000 assays	G7892
CytoTox 96® Non-Radioactive Cytotoxicity Assay (LDH, colorimetric)*	1,000 assays	G1780
GSH-Glo™ Glutathione Assay	10ml	V6911
	50ml	V6912
GSH/GSSG-Glo™ Assay	10ml	V6611
	50ml	V6612

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产品	规格	目录号
GloMax®-Multi+ Detection System with Instinct™ Software: Base Instrument with Shaking	1 each	E8032
GloMax®-Multi+ Detection System with Instinct™ Software: Base Instrument with Heating and Shaking	1 each	E9032
GloMax®-Multi+ Luminescence Module	1 each	E8041

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产品	规格	目录号
Caspase-Glo [®] 2 Assay*	10ml	G0940
Caspase-Glo [®] 6 Assay*	10ml	G0970
Caspase-Glo [®] 3/7 Assay*	2.5ml	G8090
Caspase-Glo [®] 8 Assay*	2.5ml	G8200
Caspase-Glo [®] 9 Assay*	2.5ml	G8210
Apo-ONE [®] Homogeneous Caspase-3/7 Assay	1ml	G7792
DeadEnd [™] Fluorometric TUNEL System	60 reactions	G3250
DeadEnd [™] Colorimetric TUNEL System	20 reactions	G7360
Anti-ACTIVE [®] Caspase-3 pAb	50µl	G7481
Anti-PARP p85 Fragment pAb	50µl	G7341
Anti-pS473 Akt pAb	40µl	G7441
Caspase Inhibitor Z-VAD-FMK, 20mM	50µl	G7231
	125µl	G7232

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^(a)U.S. Pat. Nos. 6,602,677 and 7,241,584, European Pat. No. 1131441, Japanese Pat. Nos. 4537573 and 4520084 and other patents pending.

^(b)U.S. Pat. No. 7,741,067, Japanese Pat. No. 4485470 and other patents pending.

^(c)U.S. Pat. No. 7,700,310, European Pat. No. 1546374 and other patents pending.

^(d)U.S. Pat. Nos 7,083,911, 7,452,663 and 7,732,128, European Pat. No. 1383914 and Japanese Pat. Nos. 4125600 and 4275715.

© 2001–2015 Promega Corporation. All Rights Reserved.

Anti-ACTIVE, Apo-ONE, Caspase-Glo, CellTiter 96, CellTiter-Blue, CellTiter-Glo, CytoTox 96 and GloMax are registered trademarks of Promega Corporation. ApoTox-Glo, ApoLive-Glo, CellTiter-Fluor, CytoTox-Fluor, CytoTox-Glo, CytoTox-ONE, DeadEnd, GSH-Glo, GSH/GSSG-Glo, Instinct, ToxGlo and Ultra-Glo are trademarks of Promega Corporation.

alamarBlue is a registered trademark of Trek Diagnostic Systems, Inc.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.