Determining Genome Targeting Efficiency using T7 Endonuclease I (M0302)

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Authenticase[®] (NEB #M0689) now available for improved mismatch detection assay performance.

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T7 Endonuclease I is also available in a convenient kit format: **EnGen Mutation Detection Kit (NEB #E3321)**

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Protocols.io also provides an **interactive version of this protocol** where you can discover and share optimizations with the research community.

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Overview: 概述：

T7 Endonuclease I recognizes and cleaves non-perfectly matched DNA. This protocol describes how to determine genome targeting efficiency by digesting annealed PCR products with T7 Endonuclease I. In the first step PCR products are produced from the genomic DNA of cells whose genomes were targeted using Cas9, TALEN, ZFN etc. In the second step, the PCR products are annealed and digested with T7 Endonuclease I. Fragments are analyzed to determine the efficiency of genome targeting.

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- a. gDNA from targeted cells (e.g. Cas9, or TALEN transfected cells)

gDNA from targeted cells （例如 Cas9 或 TALEN 转染细胞）
- b. gDNA from negative control cells (e.g. non-specific DNA transfected cells)

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PCR using Q5 High-Fidelity DNA Polymerase

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2. Gently mix the reaction. Collect all liquid to the bottom of the tube by a quick spin if necessary. Transfer PCR tubes to a PCR machine and begin thermocycling.
- 轻柔搅拌反应。如有必要，通过快速离心将所有液体收集到管底。将 PCR 管转移到 PCR 仪上并开始热循环。

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*Use of the NEB Tm Calculator is highly recommended.

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Note: Q5 Hot Start High-Fidelity 2X Master Mix does not require a separate activation step. Standard Q5 cycling conditions are recommended.

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T7 Endonuclease I digestion:

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1. Assemble reactions as follows
- 按照以下步骤组装反应

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2. Anneal the PCR products in a thermocycler using the following conditions:

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3. Add T7 Endonuclease I to the annealed PCR products

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4. Stop the reaction by adding 1.5 µl of 0.25 M EDTA.

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5. Purify the reaction using 36 µl of Ampure XP beads according to the manufacturer’s suggestion. This step is optional since 1 µl of the reaction will not interfere with analysis on an Agilent Bioanalyzer using DNA1000 reagents.

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6. Elute the DNA fragments in 20 µl of water, recovering 15 µl.

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Analysis: 分析：

- Analyze the fragmented PCR products and determine the percent of nuclease-specific cleavage products (fraction cleaved)

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- Calculate the estimated gene modification using the following formula:

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% gene modification = 100 x (1 – (1- fraction cleaved)^{1/2})

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Guschin, D.Y., et. al.(2010) A rapid and general assay for monitoring endogenous gene modification. *Methods Mol Biol*, 649, 247–256.

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Links to this resource

Product Categories: HF®, Nicking, Master Mix, Time-Saver & Other Products, Other Products, 基因编辑产品, | More +

Applications: Genome Editing Applications, Cloning & Synthetic Biology

Related Products: T7 核酸内切酶 I